



INTRAVASCULAR THERAPIES  
Ultrasound system



**vysion XS**

The ultraportable HD  
ultrasound unit

*Designed for Vascular Access*



Value Life

## vygion XS: Small and powerful for vascular access

### Perfectly suited for central venous catheter insertion

- Linear probe from 5 to 12 MHz for a visualization up to 7.4cm deep

### High-quality image

- Increased procedural accuracy
- Advanced imaging features for exceptional image quality

### Intuitive design, simplified settings

- 7" HD touch screen
- Maximum automation (automatic focalisation and frequency according to the chosen depth to obtain the best image quality)
- Ideal for real time central venous catheter placement

### Rapid boot-up time


- Start in less than 10 seconds

### Light and portable

- Only 800g (probe included)
- Autonomy up to 1h30 in continuous scanning mode




vygion XS offers the following imaging modes:




### BW

Dynamic representation of tissue



### Color

Vascularisation mapping




### Pulse wave doppler

Direction and speed of blood flow

## Key features

### Presets dedicated





### Optimized visualization

- 0.9 cm to 2.3 cm
- 2.1 cm to 5 cm
- 1.5 cm to 3.5 cm
- 3.2 cm to 7.4 cm

### Caliper

#### Accurate catheter choice

Measurement of the vein diameter and suggestion of catheter size



### Color median line

#### Practical & intuitive

Each color match the color of the needle guide associated



## Ultrasound unit: **vysion XS**

Product code	Description
<b>U-LITE</b>	Ultraportable ultrasound with a linear probe connected and fixed
<b>TABLESTAND</b>	<b>vysion XS</b> support

## Accessories

Colored median line solution works in association with:









**Reusable**  
**Vascular Access Bracket**



**Needle guide**

- 6 Depths from 0.5 cm to 3.5 cm
- For 18G to 21 G puncture needle
- 10 units/box




	Product code	Description
	<b>6261SNR</b>	<b>Reusable Vascular Access Bracket</b>
	<b>5404SNR</b>	<b>Vascular Access Bracket Needle Guide - 0.5 cm</b>
	<b>5403SNR</b>	<b>Vascular Access Bracket Needle Guide - 1.0 cm</b>
	<b>5402SNR</b>	<b>Vascular Access Bracket Needle Guide - 1.5 cm</b>
	<b>5405SNR</b>	<b>Vascular Access Bracket Needle Guide - 2.0 cm</b>
	<b>5406SNR</b>	<b>Vascular Access Bracket Needle Guide - 2.5 cm</b>
	<b>5401SNR</b>	<b>Vascular Access Bracket Needle Guide - 3.5 cm</b>

## INTRAVASCULAR THERAPIES

**For further information, please contact: [marketingbenelux@vygon.com](mailto:marketingbenelux@vygon.com)**

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

SA Vygon NV • Haachtsesteenweg 1650 Chaussée de Haecht • 1130 BRUSSELS • BELGIUM  
 Tel: +32 (0)2 706 09 50 • [www.vygon.be](http://www.vygon.be)  
 Vygon Nederland B.V. • Kerkhofstraat 21 • 5554 HG VALKENSWAARD • NEDERLAND  
 Tel: +31 (0)40 208 93 80 • [www.vygon.nl](http://www.vygon.nl)  
 VYI9117 •  Vygon Benelux

