



VASCULAR ACCESS

Haemodynamic monitoring



most-care^{Up}

Up to the beat
Product Overview



Value Life

MostCare^{Up} - Monitoring

Code	Description
VMB08MC0202E0V	MostCareUp MONITOR - On Demand version with rigid case
VMB08MC0255E0V	MostCareUp MONITOR - Endless version with rigid case








MostCare^{Up} - Cards

Code	Description
VMB04MCU1D02	MostCareUp CARD - 24 hours (1 day)
VMB04MCU1M02	MostCareUp CARD - 30 days (1 month)
VMB04MCU1Y02	MostCareUp CARD - 365 days (1 year)
VMB04MCU01U02	MostCareUp CARD - 1 use (72h)
VMB04MCU10U02	MostCareUp CARD - 10 uses (72h each)
VMB04MCU50U02	MostCareUp CARD - 50 uses (72h each)

MostCare^{Up} - Accessories

Code	Description
VMARS002505	MostCareUp TROLLEY with basket

MostCare^{Up} - Cables

Direct to transducer cable code*	Code cable Y**	Photo cable	Transducer type connection**
VMC03MUBBR	VMC03MUYBBR		BBraun transducer
VMC03MUBD	VMC03MUYBD		BD transducer
VMC03MUBIO	VMC03MUYBIO		Biosensor Utah transducer
VMC03MUDPT	VMC03MUYDPT		Codan DPT transducer
VMC03MUEDW	VMC03MUYEDW		Edwards transducer
VMC03MUMED	Not available		Medex LogiCal transducer
VMC03MUMED2	VMC03MUYMED2		Medex TranSt. transducer


* MostCareUp BP cable: direct connection to the arterial blood pressure transducer (different brand, different cables)

** MostCareUp BP Y-cable: simultaneous connection to the arterial blood pressure transducer (different brand, different cables) and the main anesthesia monitor in use. This cable allows to use only one BP transducer and arterial line to read the same blood pressure signal both in MostCareUp and the standard bedside monitor.

CRITICAL CARE

For further information, please contact : marketingbenelux@vygon.com

The specifications given in this brochure are for information only and are not, under any circumstances, of a contractual nature.

SA Vygon NV • Haachtsesteenweg 1650 Chaussée de Haecht • 1130 BRUSSELS • BELGIUM
 Tel: +32 (0)2 706 09 50 • www.vygon.be
 Vygon Nederland B.V. • Kerkhofstraat 21 • 5554 HG VALKENSWAARD • NEDERLAND
 Tel: +31 (0)40 208 93 80 • www.vygon.nl
 VY18068 •  Vygon Benelux