



ONCOLOGY
Catalogue



Qimono
by VYGON


The **high security** connection



Value Life

Preparation devices for medications/solutions

Vial

Brand name	Code	Product	Description	Indication
qimocap Ø14	7219.01 7219.010		qimocap is a vial access cap fitted with qimo ♀. It has been specifically designed to create a permanent connection with the medication vial.	Withdrawal and reconstitution of small volumes of medication in a completely closed environment.
qimocap Ø20	7219.02 7219.020			

Syringe/Extension tube

Brand name	Code	Product	Description	Indication
qimo ♂	721w		qimo ♂ is a disinfectable closed needleless safety connector. It is a safety captive adaptor for male Luer-lock hubs. Once attached, it cannot be disconnected. It connects only to qimo ♀.	To secure all the male Luer-lock fittings.
qimotransfer	7220.01		Female Luer-lock to female Luer-lock connector, with qimo ♀ on both ends of the connector.	To transfer medications from syringe to syringe.

Priming device

Brand name	Code	Product	Description	Indication
qimoprime	7210.81		qimoprime is a priming device which can be connected to any product fitted with qimo ♂.	Priming device for qimo ♂

Elastomeric Pump

Infusion line

Brand name	Code	Product	Description	Indication
qimoprime & qimo ♂	7210.01		Set supplied with qimo ♂ and qimoprime in the same packaging.	Can be used to secure distal end of elastomeric pump tubing. Need to be connected to a qimo ♀

Patient's vascular access

Brand name	Code	Product	Description	Indication
qimo ♀	7210.02		qimo ♀ is a multi-purpose, closed, needleless and disinfectable safety connector, supplied in individual packaging.	qimo ♀ can be used to connect the tubing of the distal end of the elastomeric pump secured with qimo ♂, to access an implantable port or any other patient vascular device.

Patient's vascular access & Infusion line

Brand name	Code	Product	Description	Indication
qimoprime, qimo ♂ & qimo ♀	7210.04		Set with qimoprime, qimo ♂ and qimo ♀ in the same packaging.	Can be used to secure the filling port of the elastomeric pump and the distal end of its tubing, but other applications are possible.

● Infusion devices


Infusion line

Brand name	Code	Product	Description	Indication
qimoprime & qimo ♂	7210.01		Set supplied with qimo ♂ and qimoprime in the same packaging.	To secure both, the male Luer-lock and the female Luer-lock end of infusion devices.

Patient's vascular access

Brand name	Code	Product	Description	Indication
qimo ♀	7210.02		qimo ♀ is a multi-purpose, closed, needleless and disinfectable safety connector.	To connect to any other vascular access, Huber needles, CVC... safely and with a closed system

● Vesical instillation

Brand name	Code	Product	Description	Indication
qimoadapt	7210.06		qimoadapt is a safe adaptor for the connection of a syringe on the funnel of a tube. qimo ♀ allows to connect a syringe secured with qimo ♂.	Safe bladder instillation

INTRAVASCULAR THERAPIES

For further information, please contact: marketingbenelux@vygon.com
The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

SA Vygon NV • Haachtsesteenweg 1650 Chaussée de Haecht • 1130 BRUSSELS • BELGIUM
Tel: +32 (0)2 706 09 50 • www.vygon.be
Vygon Nederland B.V. • Kerkhofstraat 21 • 5554 HG VALKENSWAARD • NEDERLAND
Tel: +31 (0)40 208 93 80 • www.vygon.nl
VY19101 • [in](#) Vygon Benelux



DB ONCO 18 416 E - 07/19