



ENTERAL FEEDING CATALOGUE



Want to know more about safe
feeding system ?


www.safe-enteral.com

www.vygon.com

Vygon in the world

OUR STORY

Pierre Simonet founded Vygon in France in 1962. He was one of the first to introduce the concept of disposable medical equipment. Shortly thereafter, the company established subsidiaries in Europe, the United States, India, South America, and more. Today it is present worldwide. Vygon also relies on many production units, starting with its first plant in France and continuing on to its latest, brand-new plant that opened in Portugal in 2015.



Vygon,
committed to
its customers
worldwide

Vygon is a company that designs, manufactures and markets a wide range of medical devices for 4000 customers worldwide. Its success is due to its teams, who are available and committed. They are at the service of caregivers to offer them products and solutions that improve the safety and performance of care. The company that invented the single-use medical device, a true revolution in hospital care practices, has built its offer on the complementary skills that make up its business: industrial and technological know-how, clinical knowledge and understanding of hospital structures, and its own resources in research and development. Today, Vygon is an international company but it also remains a family business. The stability and commitment of its shareholders allow for a sustainable development that brings projects to fruition and its value system is a guide.

VYGON OFFERS A SAFE AND SECURE ENTERAL FEEDING SYSTEM FOR EVERY PATIENT

In the 2000's, the International Organisation for Standardisation (ISO) decided to mitigate the likelihood of tubing misconnections and create the ISO 80369 series in order to increase patient safety by specifying designs of non-interconnectable connectors for liquids and gases used in healthcare applications (respiratory, enteral, urology, limb cuff inflation, neuraxial, and I.V.).

In order to get the technical and market expertise, ISO encouraged industry to be part of the working groups. That is why Vygon, with 20 years of experience in safety enteral feeding systems, joined the ISO 80369 Committee from the beginning.

In July 2016, the International Organisation for Standardisation published the ISO 80369-3 standard, which defines a standardised safety enteral feeding connection. As with **nutrisafe2**, this new enteral connection, named ENFit™, prevents the risk of misconnections (e.g. inject milk in the intravenous line) by its specific non-interconnectable design.

Patient safety, along with high quality care, is at the top of Vygon's priorities for all populations. That is why Vygon goes further than the ISO 80369 standard by proposing not only an ENFit™ range, but also a dedicated safety enteral range for neonates and newborns.



nutriFit®, A COMPREHENSIVE ENFit™ PRODUCT RANGE

- **Compliant with ISO 80369-3:** nutriFit® incorporates the standardised safety enteral connector.
- **Specific:** only dedicated to enteral application, nutriFit® prevents misconnection with other clinical applications (respiratory, I.V., etc.).
- **Secure:** with its twist design, nutriFit® avoids accidental disconnections.
- **Universal:** nutriFit® can be connected to any ENFit™ system worldwide.

Due to the large dimensions of the ENFit™ connectors, nutriFit® products are especially well adapted for the old children and adults.

nutrisafe2, A SAFETY ENTERAL PRODUCT RANGE DEDICATED TO NEONATES AND NEWBORNS

Due to the neonatal patient's tiny size, they require special care with dedicated medical devices. Therefore, Vygon decided to maintain nutrisafe2 for them:

- **Safe:** Compliant with ISO 80369-1, nutrisafe2 incorporates a specific safety enteral connector, which was assessed to prevent misconnection risks with other clinical applications (respiratory, I.V., etc.) in line with the ISO 80369 general requirements.
- **Small:** Specially designed for preterm infants, nutrisafe2 compact design minimises the size, weight and dead space of the connection.
- **Accurate:** Thanks to its compact design, nutrisafe2 ensures high accuracy in drug administration in very low volumes and reduces the overdosing risk by a factor of 4 when compared to ENFit™ or Low Dose Tip (LDT) connections.
- **More than 10 years of experience in neonatology:** Vygon took on board customer feedback to continuously improve its nutrisafe2 range for neonatal use and validated its safety aspect.



TABLE OF CONTENTS

Pumps And Sets

page 8

nutrisafe2
The safe enteral
feeding system for
neonates

page 19

| | |
|--|----|
| ENTERAL SYRINGE PUMP | 9 |
| ENTERAL FEEDING PUMPS | 10 |
| easy lit | 10 |
| easy moov6 | 11 |
| Giving sets for pumps | 12 |
| easy lit range giving pump sets | 13 |
| easy moov6 giving pump sets | 14 |
| Magnetic giving pump sets | 15 |
| GRAVITY GIVING SETS | 16 |
| PUMP ACCESSORIES | 17 |
| Bags for pump | 17 |
| Upper double line adaptors | 17 |
| Universal bag spike/clamp | 18 |
| ENPlus bag adaptors | 18 |
| easy moov6 accessories | 18 |
| FEEDING TUBES | 20 |
| PVC feeding tubes | 20 |
| Polyurethane feeding tubes | 21 |
| Silicone feeding tubes | 24 |
| GASTRIC ASPIRATION | 26 |
| Dual flow gastric tube | 26 |
| SYRINGES | 27 |
| Transparent syringes | 27 |
| Amber coloured syringes | 28 |
| Re-usable syringes | 28 |
| Syringe trays | 29 |
| ACCESSORIES | 30 |
| Extension tubes | 30 |
| Syringe accessories | 32 |
| Sampling devices | 33 |
| Other accessories | 36 |

ENFit™ / nutrifit®

Safe enteral feeding
systems for children
and adults

page 39

Non-safety enteral feeding devices

page 55

ENFit™ / nutrifit®

Transition period
accessories

page 67

| | |
|---|----|
| FEEDING TUBES | 40 |
| PVC feeding tubes | 40 |
| Polyurethane feeding tubes | 42 |
| Silicone feeding tubes | 46 |
| SYRINGES | 48 |
| Low Dose Tip syringes | 48 |
| Transparent syringes | 49 |
| Re-usable syringes | 49 |
| SYRINGE TRAYS | 50 |
| Accessories | 51 |
| Extension tubes | 51 |
| Syringes accessories | 52 |
| Sampling devices | 53 |
| Other accessories | 54 |
| FUNNELED FEEDING TUBES | 57 |
| PVC feeding tubes | 57 |
| Polyurethane feeding tubes | 59 |
| Silicone feeding tubes | 60 |
| JEJUNOSTOMY TRAY | 62 |
| GASTRIC ASPIRATION | 63 |
| Faucher gastric tubes | 63 |
| Dual flow gastric tubes | 64 |
| ACCESSORIES | 66 |
| Stylet for silicone feeding tubes | 66 |
| Non-return valve for dual flow gastric tube | 66 |
| CHRISTMAS TREE ADAPTORS | 68 |
| FUNNEL ADAPTORS | 68 |

PUMPS AND SETS

ENTERAL SYRINGE PUMP

Designed to meet all patients needs in stationary medical care, this easy-to-use high-quality syringe pump is dedicated to accurate enteral feeding.

An indicator clearly indicates the administration status. As soon as a syringe is placed, its volume is automatically identified (10, 20, 30/35 or 50ml). The syringe calibration is compatible with A-VY, C-GON, VY-P and **nutrifit** syringes (see table below). Buttons have symbols to allow a better identification.

The syringe pump is fixed with a multi-direction clamp that enables a vertical or horizontal pole installation. The device can be stabilized thanks to an anti-slip mat present at the bottom.

- Features & benefits:**
- Compatible with all brands of **nutrisafe2** and **nutrifit** syringes
 - Accuracy of flowrate including syringe $\pm 2\%$
 - Optic and sound alarms
 - Possibility to adjust sound level
 - Up to 2000 history events
 - Delivered volume inquiry
 - Occlusion detection pressure setting
- Weight: 2.25kg (with clamp base)
Purge rate: 400 to 1500ml/h
Flow rate: 0.1 to 1500 ml/h
Battery: Ni-MH
Battery autonomy: 6 hours



The following syringes are already calibrated on the device. Others channels are available to calibrate new syringe brands (contact your biomedical department).

| | nutrisafe2 | | | nutrifit |
|--------------|-------------------|------------------|------------------|-----------------------------|
| Syringe vol. | Channel 2 / C-GON | Channel 3 / A-VY | Channel 4 / VY-P | Channel 5 / nutrifit |
| 10ml | 1015.102 | 1015.103 | 1015.104 | 1015.107 |
| 20ml | 1015.202 | 1015.213 | 1015.204 | 1015.207 |
| 30 / 35ml | 1015.352 | - | - | - |
| 50ml | 1015.602 | 1015.603 | 1015.604 | 1015.607 |

- Supplied with:
- Power cords Europe (plug C) / UK (plug G) / US (plug A)
 - User manual and sticker
 - Packing list
 - Certificate of compliance
 - Warranty card

0VKSN50C6
Syringe pump

ENTERAL FEEDING PUMPS

easy**lit**6

easy**lit**6 is an enteral feeding pump, designed for stationary use, ideal for overnight feeding, hospitalised patients and elderly people.

easy**lit**6 is a stationary pump allowing flowrate and volume programming.

Features & benefits:

- Easy to use: short learning curve
- Powerful: the rotary peristaltic motor helps delivering high viscosity food without any obstruction
- Quiet: the continuous rotary peristaltic system avoids repetitive start and stop noises
- Safe: visual and audible alarms, calibrated to avoid untimely alarms due to overpressure (in case of cough, brief kinkings, etc.).



- Weight: 1400gr (3.08lb)
- Flow rate: 1 - 600ml/h in 1ml increment
- Volume: 1 to 9999ml in 1ml increment
- Battery: Ni-MH
- Battery autonomy: 7h at 125ml/h
- Selfpriming

OVEPL6

easy**lit**6 stationary pump

easy**moov**6

easy**moov**6 is an enteral feeding pump, designed to meet all ambulatory patient needs thanks to its innovative design and ergonomic shape.

The new easy**moov**6 ambulatory pump is also suitable for stationary use thanks to its 180 degree pole clamp. This pump offers comfort for patients and medical professionals to meet all their needs in hospital and homecare.

Intuitive and available in 8 languages, this new pump has a large screen showing useful and relevant information during feeding such as flow rate, volume, time left before the end of the feeding and total volume already administered.

A Menu option is also available, allowing patients, caregivers or medical professionals to change the pump settings, to set two different nutrition modes and to see the pump data without any software.

Features & benefits:

- Intuitive and easy to use: short learning curve. All the information is available in a large screen, allowing the user to manage it quickly
- Powerful: the rotary peristaltic motor helps delivering high viscosity food without any obstruction, guaranteeing silent and constant flow rate throughout the feeding
- Quiet: the continuous rotary peristaltic system avoids repetitive start and stop noises allowing patient feeding during the night without being disturbed
- Safe: visual and audible alarms, calibrated to avoid untimely alarms due to overpressure (in case of cough, brief kinkings, etc.).



- Weight: 650gr (with battery)
- Flow rate: 1 - 400ml/h in 1ml increment
- Volume: 1 to 9999ml in 1ml increment
- Battery: Ni-MH
- Battery autonomy: 20h at 125ml/h
- Selfpriming

OVEPM6



easy**moov**6 ambulatory pump

GIVING SETS FOR PUMPS

A full range of giving sets designed to cover all healthcare professionals' requirements and to improve patient safety and comfort.



The range includes giving sets with or without bag, single or double lines and double line with one line dedicated to hydration and fitted with a non-return valve.

The range is designed with the safest components available ensuring the safest delivery of nutrition:






Proximal end

Fitted with the enteral safety bag connector compatible with ENPlus bags to prevent any risk of misconnection with parenteral nutrition bags.



Lateral port

Fitted with the safe ENFit™ or nutrisafe2 connector compatible with ENFit™ or nutrisafe2 syringes.



Distal end

Includes an ENFit™ or nutrisafe2 connector with a removable christmas tree adaptor.

easylit RANGE GIVING PUMP SETS

This range of giving pump sets is compatible with easylit pumps and **also with P7000, Compat, Flocare 800 pumps**. Giving pump sets fitted with a non-return valve are dedicated to hydration.



| easylit giving sets | | | | | | |
|----------------------------|----------------|------------------------|--------------------|---------------------------------------|------------------------|---------|
| | Code | Proximal connection | Lateral connection | Distal connection | Additional information | Qty/box |
| Pump sets with one line | 0VENEJ34S | EN+ | ENFit™ | ENFit™ & Christmas tree connector | - | 30 |
| | 0VENEJ34SCS | Screw cap | | | - | 30 |
| | 0VENP2GBJ34S | 1200ml bag | | | - | 50 |
| | 0VENEJ34V | EN+ | nutrisafe2 | nutrisafe2 & Christmas tree connector | - | 30 |
| | 0VENEJ34VCS | Screw cap | | | - | 30 |
| | 0VENP1GBJ34V | 500ml bag | | | - | 50 |
| | 0VENP2GBJ34V | 1200ml bag | | | - | 50 |
| Pump sets with double line | 0VENEEJ24S | EN+ ; EN+ | ENFit™ | ENFit™ & Christmas tree connector | - | 30 |
| | 0VENEEJ24SH | EN+ ; EN+ | | | With non return valve | 30 |
| | 0VENEP2J24SHCS | 1200ml bag ; screw cap | | | With non return valve | 50 |
| | 0VENEEJ24V | EN+ ; EN+ | nutrisafe2 | nutrisafe2 & Christmas tree connector | - | 30 |

For more giving pump sets, please contact your Vygon representative.

easy**moov**6 GIVING PUMP SETS

This range of giving pump sets is compatible with easy**moov**6 only. Giving pump sets fitted with a non-return valve are dedicated to hydration.



| easy moov 6 giving pump sets | | | | | |
|-------------------------------------|-----------------------|--------------------|---------------------------------------|---------------------------|----------|
| Code | Proximal connection | Lateral connection | Distal connection | Additional information | Qty/ box |
| 0VENEP14S | EN+ | ENFit™ | ENFit™ & Christmas tree connector | - | 30 |
| 0VENEP08S | EN+ | | | With non return valve | 30 |
| 0VENEP14SCS | Screw cap | | | - | 30 |
| 0VENFP14S | Spike | | | - | 30 |
| 0VENEP14R | EN+ | | | With 3 way stopcock | 30 |
| 0VENEP14RCS | Screw cap | | | With 3 way stopcock | 30 |
| 0VENP1P14S | 500ml bag | | | - | 50 |
| 0VENP2P14S | 1200ml bag | | | - | 50 |
| 0VENEP11S | EN+ | - | ENFit™ | - | 30 |
| 0VENEP11SCS | Screw cap | - | ENFit™ | - | 30 |
| 0VENEP14V | EN+ | nutrisafe2 | nutrisafe2 & Christmas tree connector | | 30 |
| 0VENEP08V | EN+ | | | With non return valve | 30 |
| 0VENEP14VCS | Screw cap | | | - | 30 |
| 0VENP1P14V | 500ml bag | | | - | 50 |
| 0VENEFP14S | EN+ ; EN+ | ENFit™ | ENFit™ & Christmas tree connector | - | 30 |
| 0VENEFP14SH | EN+ ; EN+ | | | With anti free-flow valve | 30 |
| 0VENFEP14S | Spike ; EN+ | | | - | 30 |
| 0VENEP2P14SHCS | 1200ml bag; screw cap | | | With anti free-flow valve | 50 |
| 0VEN2P2P14SH | 2 x 1200ml bags | | | With anti free-flow valve | 50 |

For more giving pump sets, please contact your Vygon representative.

MAGNETIC GIVING PUMP SETS

This range of giving pump sets is compatible with magnetic detection pumps such as K224, K324, K624 and K2100. Giving pump sets fitted with a non-return valve are dedicated to hydration.



| Magnetic giving pump sets | | | | | | |
|----------------------------|--------------|---------------------|--------------------|-----------------------------|------------------------|----------|
| | Codes | Proximal connection | Lateral connection | Distal connection | Additional information | Qty/ box |
| Pump sets with one line | 0VENEK34S | EN+ | ENFit™ | ENFit™ & Christmas tree | - | 30 |
| | 0VENEK34SCS | Screw cap | | | - | 30 |
| | 0VENP2GBK34S | 1200ml bag | | | - | 50 |
| | 0VENEK34V | EN+ | nutrisafe2 | nutrisafe2 & Christmas tree | - | 30 |
| Pump sets with double line | 0VENEK24S | EN+ ; EN+ | ENFit™ | ENFit™ & Christmas tree | - | 30 |
| | 0VENEK24SH | EN+ ; EN+ | | | With non return valve | 30 |

For more giving pump sets, please contact your Vygon representative.

GRAVITY GIVING SETS

This range of giving sets is designed for gravity use only. It includes two different diameters: a small one and a large one, depending on viscosity and flow rate.



| Gravity giving sets | | | | | |
|---------------------|------------|---------------------|--------------------|-------------------------|---------|
| | Codes | Proximal connection | Lateral connection | Distal connection | Qty/box |
| Small diameter | 0VENE74S | EN+ | ENFit™ | ENFit™ & Christmas tree | 30 |
| | 0VENE74SCS | Screw cap | | | 30 |
| | 0VENE74GCS | Screw cap | | | 30 |
| Large diameter | 0VENE94S | EN+ | | | 30 |
| | 0VENE94SCS | Screw cap | | | 30 |

For more gravity giving sets, please contact your Vygon representative.

PUMP ACCESSORIES


BAGS FOR PUMP

- All bags feature:
- The new enteral safety bag port compatible ENPlus
 - An ice compartment for cooling the enteral nutrition
 - Possibility to indicate patient name and information about nutrition
 - A large orifice to pour liquids & easy to close
 - Non progressive Roberts clamp

| Code | Designation | Quantity/box |
|----------|--------------------|--------------|
| 0VENP1GT | 500ml enteral bag | 50 |
| 0VENP2GT | 1200ml enteral bag | 50 |
| 0VENP3GT | 1500ml enteral bag | 50 |



UPPER DOUBLE LINE ADAPTORS



| Code | Designation | Quantity/box |
|-----------|---|--------------|
| 0VEACC127 | Upper double line adaptor fitted with ENPlus connection | 100 |
| 0VEACC126 | Upper double line adaptor fitted with ENPlus connection ; with hydration line | 100 |

UNIVERSAL BAG SPIKE/CLAMP

| Code | Designation | Quantity/box |
|-----------|-------------------------------|--------------|
| 0VEACC015 | Bag spike | 300 |
| 0VEACC023 | Removable clamp (non-sterile) | 50 |



ENPLUS BAG ADAPTORS



| Code | Designation | Quantity/box |
|--------------|-------------------------------------|--------------|
| 0VEACC114 | ENPlus bag adaptor (non-sterile) | 100 |
| 0VTACC300ENT | Screw cap | 100 |
| 0VEACC360 | ENPlus / ENFit connector with clamp | 100 |

easy**moov**6 ACCESSORIES



| Code | Designation | Quantity/box |
|-------------|---|--------------|
| 0VEPM6BAG | easy moov 6 compact back pack | 1 |
| 1MWPM6CHAx | easy moov 6 power charger | 1 |
| 1MWPM6ACP | easy moov 6 pole clamp | 1 |
| 0VESUPPOMP | easy moov 6 pump table stand | 1 |
| 1MWEASYSOFT | easy soft for easy moov 6 | 1 |

x: depending on the type of plug required

nutrisafe2

THE SAFE ENTERAL FEEDING SYSTEM
FOR NEONATES

FEEDING TUBES

PVC FEEDING TUBES

Especially designed for neonatal patients and newborns, the PVC nutriSafe2 feeding tubes are intended to be used for short-term (up to 5 days) and feature the following characteristics:

- Made of PVC
- Closed distal tip with 2 lateral eyes (4 to 12Fr) and 4 lateral eyes (14 to 18Fr)
- Small cap with colour code to quickly identify the size of the tube
- Numerical markings every centimetre:
 - from 5cm to 35cm for codes 364, 361 & 363
 - from 5cm to 70cm for code 362
- Radiopaque
- Transparent to better visualize the administered solution.

| Code | Tube | | | | | Colour code | Quantity |
|---------|-----------|----|-----------|-----------|---------------|-------------|-----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | | Box/ case |
| 364.042 | 40 | 4 | 1.5 | 0.8 | 0.24 | Purple | 25/600 |
| 364.052 | 40 | 5 | 1.7 | 1.0 | 0.36 | Grey | 25/600 |
| 364.062 | 40 | 6 | 2.0 | 1.2 | 0.51 | Green | 25/600 |
| 364.082 | 40 | 8 | 2.5 | 1.5 | 0.78 | Blue | 25/600 |
| 361.042 | 50 | 4 | 1.5 | 0.8 | 0.28 | Purple | 50/600 |
| 361.052 | 50 | 5 | 1.7 | 1.0 | 0.43 | Grey | 50/600 |
| 361.062 | 50 | 6 | 2.0 | 1.2 | 0.60 | Green | 50/600 |
| 361.082 | 50 | 8 | 2.5 | 1.5 | 0.93 | Blue | 50/600 |
| 361.102 | 50 | 10 | 3.3 | 2 | 1.62 | Black | 50/600 |
| 363.042 | 75 | 4 | 1.5 | 0.8 | 0.41 | Purple | 25/600 |
| 363.052 | 75 | 5 | 1.7 | 1.0 | 0.62 | Grey | 25/600 |
| 363.062 | 75 | 6 | 2.0 | 1.2 | 0.89 | Green | 25/600 |
| 363.082 | 75 | 8 | 2.5 | 1.5 | 1.37 | Blue | 25/600 |
| 363.102 | 75 | 10 | 3.3 | 2 | 2.41 | Black | 25/600 |
| 363.122 | 75 | 12 | 4 | 2.5 | 3.75 | White | 25/600 |
| 362.042 | 125 | 4 | 1.5 | 0.8 | 0.66 | Purple | 25/600 |
| 362.052 | 125 | 5 | 1.7 | 1.0 | 1.02 | Grey | 25/600 |
| 362.062 | 125 | 6 | 2.0 | 1.2 | 1.45 | Green | 25/600 |
| 362.082 | 125 | 8 | 2.5 | 1.5 | 2.25 | Blue | 25/600 |
| 362.102 | 125 | 10 | 3.3 | 2.0 | 3.98 | Black | 25/600 |
| 362.122 | 125 | 12 | 4 | 2.5 | 6.20 | White | 25/300 |
| 362.142 | 125 | 14 | 4.5 | 3 | 8.56 | Dark green | 50/200 |
| 362.162 | 125 | 16 | 5 | 3.5 | 11.96 | Orange | 50/200 |
| 362.182 | 125 | 18 | 6 | 4 | 16.12 | Red | 50/200 |

Feeding tubes are also available in orange.



POLYURETHANE FEEDING TUBES

WITHOUT STYLET

Designed for neonatal patients and newborns, the PUR nutriSafe2 feeding tubes are intended to be used for medium-term (1 to 4 weeks) and feature the following characteristics:

- Made of PUR
- Closed distal tip with 2 lateral eyes (4 to 12Fr) and 4 lateral eyes (14 to 18Fr)
- Small cap with colour code to quickly identify the size
- Numerical markings every centimetre:
 - from 5cm to 35cm for codes 1364, 1361 & 1363
 - from 5cm to 70cm for codes 1362 & 1365
- Radiopaque line from 4Fr to 6Fr, and ¾ radiopaque tubing from 8Fr
- Transparent to better visualize the administered solution.




| Code | Tube | | | | | Colour code | Quantity |
|----------|-----------|----|-----------|-----------|---------------|-------------|-----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | | Box/ case |
| 1364.042 | 40 | 4 | 1.5 | 0.8 | 0.26 | Purple | 25/600 |
| 1364.052 | 40 | 5 | 1.7 | 1 | 0.38 | Grey | 25/600 |
| 1364.062 | 40 | 6 | 2 | 1.2 | 0.53 | Green | 25/600 |
| 1364.082 | 40 | 8 | 2.5 | 1.5 | 0.8 | Blue | 25/600 |
| 1361.042 | 50 | 4 | 1.5 | 0.8 | 0.31 | Purple | 50/600 |
| 1361.052 | 50 | 5 | 1.7 | 1 | 0.46 | Grey | 50/600 |
| 1361.062 | 50 | 6 | 2 | 1.2 | 0.64 | Green | 50/600 |
| 1361.082 | 50 | 8 | 2.5 | 1.5 | 0.97 | Blue | 50/600 |
| 1361.102 | 50 | 10 | 3.3 | 2 | 1.7 | Black | 50/600 |
| 1363.042 | 75 | 4 | 1.5 | 0.8 | 0.44 | Purple | 25/600 |
| 1363.052 | 75 | 5 | 1.7 | 1 | 0.66 | Grey | 25/600 |
| 1363.062 | 75 | 6 | 2 | 1.2 | 0.92 | Green | 25/600 |
| 1363.082 | 75 | 8 | 2.5 | 1.5 | 1.4 | Blue | 25/600 |
| 1363.102 | 75 | 10 | 3.3 | 2 | 2.5 | Black | 25/600 |
| 1363.122 | 75 | 12 | 4 | 2.5 | 3.8 | White | 25/300 |
| 1362.042 | 125 | 4 | 1.5 | 0.8 | 0.69 | Purple | 25/600 |
| 1362.052 | 125 | 5 | 1.7 | 1 | 1.05 | Grey | 25/600 |
| 1362.062 | 125 | 6 | 2 | 1.2 | 1.5 | Green | 25/600 |
| 1362.082 | 125 | 8 | 2.5 | 1.5 | 2.3 | Blue | 25/600 |
| 1362.102 | 125 | 10 | 3.3 | 2 | 4 | Black | 25/600 |
| 1362.122 | 125 | 12 | 4 | 2.5 | 6.3 | White | 25/300 |
| 1362.142 | 125 | 14 | 4.5 | 3 | 9 | Dark green | 50/200 |
| 1362.162 | 125 | 16 | 5 | 3.5 | 12.3 | Orange | 50/200 |
| 1362.182 | 125 | 18 | 6 | 4 | 16 | Red | 20/200 |
| 1365.042 | 160 | 4 | 1.5 | 0.8 | 0.86 | Purple | 25/600 |
| 1365.062 | 160 | 6 | 2 | 1.2 | 1.9 | Green | 25/600 |
| 1365.082 | 160 | 8 | 2.5 | 1.5 | 2.9 | Blue | 25/600 |

HIGH FLOW WITHOUT STYLET

Designed for neonatal patients needing a higher flowrate, the 4Fr High Flow polyurethane feeding tube enables the flow rate to be increased by +48%⁽¹⁾. It is intended to be used up to 29 days and features the following characteristics:

- Made of PUR
- Closed distal tip with 2 lateral eyes
- Small cap with colour code to quickly identify the size
- Numerical markings every centimetre from 5cm to 35cm
- Radiopaque line
- Transparent to better visualize the administered solution




| Code | Tube | | | | | Colour code | Quantity |
|-----------|-----------|----|-----------|-----------|---------------|-------------|-----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | | Box/ case |
| 1364.0425 | 40 | 4 | 1.5 | 1 | 0.38 | Purple | 25/600 |

(1)RST_CDC15004B_135_21_A1_Flow comparison of 4Fr nutri**safe**2 feeding tube

WITH STYLET

Designed for neonatal patients and newborns, the PUR nutri**safe**2 feeding tubes are intended to be used for medium-term (1 to 4 weeks) and feature the following characteristics:

- Made of PUR
- Stylet protected by a ball end to facilitate the placement
- Closed or open distal tip with 2 lateral eyes
- Small cap with colour code to quickly identify the size
- Numerical markings every centimetre:
 - from 5cm to 35cm for code 1363
 - from 5cm to 70cm for code 1362
- Radiopaque.



| Code | Tube | | | | | Distal tip | Colour code | Quantity |
|-----------|-----------|----|-----------|-----------|---------------|------------|-------------|-----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | | | Box/ case |
| 1363.0621 | 75 | 6 | 2 | 1.2 | 0.93 | Closed | Green | 25/300 |
| 1363.0821 | 75 | 8 | 2.5 | 1.7 | 1.8 | Closed | Blue | 25/300 |
| 1363.0822 | 75 | 8 | 2.5 | 1.7 | 1.8 | Open | Blue | 25/300 |
| 1362.0621 | 125 | 6 | 2 | 1.2 | 1.5 | Closed | Green | 25/300 |
| 1362.0821 | 125 | 8 | 2.5 | 1.7 | 2.9 | Closed | Blue | 25/300 |
| 1362.0822 | 125 | 8 | 2.5 | 1.7 | 2.9 | Open | Blue | 25/300 |

SILICONE FEEDING TUBES

WITHOUT WEIGHTED TIP

Designed for neonatal patients and newborns, the silicone gastro-duodenal nutri**safe**2 feeding tubes are intended to be used up to 4 weeks and feature the following characteristics:

- Made of silicone
- Closed distal tip with 2 lateral eyes
- Small cap with colour code to quickly identify the size
- Numerical markings every centimetre:
 - from 5cm to 35cm for code 2331
 - from 5cm to 70cm for code 2332
- Radiopaque
- Transparent to visualize the administered solution.



| Code | Tube | | | | | Colour code | Quantity |
|----------|-----------|----|-----------|-----------|---------------|-------------|----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | | Box/case |
| 2331.042 | 50 | 4 | 1.5 | 0.8 | 0.29 | Purple | 10/480 |
| 2331.052 | 50 | 5 | 1.7 | 1 | 0.44 | Grey | 10/480 |
| 2331.062 | 50 | 6 | 2 | 1.2 | 0.62 | Green | 10/480 |
| 2331.082 | 50 | 8 | 2.5 | 1.5 | 0.96 | Blue | 10/480 |
| 2331.102 | 50 | 10 | 3.3 | 2 | 1.68 | Black | 10/480 |
| 2332.042 | 125 | 4 | 1.5 | 0.8 | 0.66 | Purple | 10/480 |
| 2332.062 | 125 | 6 | 2 | 1.2 | 1.44 | Green | 10/480 |
| 2332.082 | 125 | 8 | 2.5 | 1.5 | 2.34 | Blue | 10/480 |
| 2332.102 | 125 | 10 | 3.3 | 2 | 4 | Black | 10/480 |

Feeding tubes are also available in orange.

SILICONE FEEDING TUBES

WITH WEIGHTED TIP

Designed for neonatal patients and newborns, the silicone gastro-duodenal nutri**safe**2 feeding tubes with weighted tip are intended to be used up to 4 weeks and feature the following characteristics:

- Made of silicone
- Weighted tip with metallic olive
- Closed distal tip with 2 lateral eyes
- Small cap with colour code to quickly identify the size
- Numerical markings every centimetre from 5cm to 70cm
- Radiopaque
- Transparent to visualize the administered solution.



| Code | Tube | | | | | Colour code | Quantity |
|----------|-----------|----|-----------|-----------|---------------|-------------|----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | | Box/case |
| 2395.062 | 125 | 6 | 2.0 | 1.0 | 1.0 | Green | 10/480 |
| 2395.082 | 125 | 8 | 2.5 | 1.3 | 2.21 | Blue | 10/480 |

GASTRIC ASPIRATION

DUAL FLOW GASTRIC TUBE

Dual flow gastric tube for continuous gastric aspiration with nutrisafe2 connection:

- Made of PVC
- The main channel is fitted with a removable taper connector (code 368.12)
- The extension tube of the secondary channel ends with a male nutrisafe2 hub
- Closed distal tip with 6 lateral eyes
- Numerical markings every centimetre from 5cm to 45cm
- Radiopaque
- The product is also supplied with a straight double taper connector
- Length: 60 cm
- Transparent to visualize the administered solution.

Available in 6Fr (code 340.062) and 8Fr (code 340.082).
Quantity box/case: 10/180



SYRINGES

TRANSPARENT SYRINGES

The transparent nutrisafe2 syringes incorporate a purple-coloured plunger and “Enteral” print, enabling instant identification for enteral use. Sterile and individually packed.



• C-GON SYRINGES

| Code | C-GON Syringes | | | Quantity |
|----------|----------------|------------------|-----------------------|----------|
| | Volume ml | Graduation every | nutrisafe2 connection | Box/case |
| 1015.002 | 0.5 | 0.01 ml | Centered, slip | 100/3200 |
| 1015.352 | 35 | 1 ml | Centered, lockable | 50/200 |



• A-VY SYRINGES

| Code | A-VY Syringes | | | Quantity |
|----------|---------------|------------------|-----------------------|----------|
| | Volume ml | Graduation every | nutrisafe2 connection | Box/case |
| 1015.013 | 1 | 0.01 ml | Centered, slip | 160/1600 |
| 1015.023 | 2.5 | 0.1 ml | Centered, slip | 130/1300 |
| 1015.053 | 5 | 0.2 ml | Centered, lockable | 100/1000 |
| 1015.103 | 10 | 0.2 ml | Offset, lockable | 100/1000 |
| 1015.213 | 20 | 0.5 ml / 1 ml | Offset, lockable | 100/600 |
| 1015.603 | 60 | 1 ml | Offset, lockable | 50/300 |

These syringes are also available in orange. Please contact your Vygon representative.



• VY-P SYRINGES

| Code | Cap | VY-P Syringes | | | Quantity |
|-----------|-----|---------------|------------------|-----------------------|----------|
| | | Volume ml | Graduation every | nutrisafe2 connection | Box/case |
| 1015.014 | No | 1 | 0.01 ml | Centered, slip | 100/2400 |
| 1015.024 | No | 2.5 | 0.1 ml | Centered, slip | 100/2000 |
| 1015.054 | No | 5 | 0.2 ml | Centered, lockable | 100/1600 |
| 1015.104 | No | 10 | 0.2 ml | Centered, lockable | 100/1000 |
| 1015.204 | No | 20 | 0.5 ml | Centered, lockable | 80/800 |
| 1015.214 | No | 20 | 0.5 ml | Offset, lockable | 80/800 |
| 1015.604 | No | 60 | 1 ml | Offset, lockable | 60/240 |
| 1015.1041 | Yes | 10 | 0.2 ml | Centered, lockable | 100/1000 |
| 1015.2141 | Yes | 20 | 0.5 ml | Offset, lockable | 80/800 |
| 1015.6041 | Yes | 60 | 1 ml | Offset, lockable | 60/240 |

AMBER COLOURED SYRINGES

The nutrisafe2 amber-coloured syringes reduce the exposure to light for medications. The syringes incorporate a purple-coloured plunger and “Enteral” print, enabling instant identification for enteral use.

These syringes are not sterile, packed by 50 units and supplied with 50 caps.



| Code | Amber-coloured syringes | | | Quantity |
|----------|-------------------------|------------------|-----------------------|----------|
| | Volume ml | Graduation every | nutrisafe2 connection | Box/case |
| 1022.012 | 1 | 0.01 ml | Centered, slip | 50/2400 |
| 1022.022 | 2.5 | 0.1 ml | Centered, slip | 50/2400 |
| 1022.052 | 5 | 0.2 ml | Centered, lockable | 50/1200 |
| 1022.102 | 10 | 1 ml | Centered, lockable | 50/450 |

RE-USABLE SYRINGES

The VY-P nutrisafe2 reusable syringes have been specifically designed for homecare enteral nutrition and medication. The syringes incorporate a purple-coloured plunger and “Enteral” print, enabling instant identification for enteral use.

These syringes are sterile (for the first use) and individually packed.



| Code | Re-usable syringes | | | Quantity |
|----------|--------------------|---------------------|-----------------------|----------|
| | Volume ml | Graduation every ml | nutrisafe2 connection | Box/case |
| 1012.014 | 1 | 0.01 | Centered, slip | 100/2400 |
| 1012.024 | 2.5 | 0.1 | Centered, slip | 100/2000 |
| 1012.054 | 5 | 0.2 | Centered, lockable | 100/1600 |
| 1012.104 | 10 | 0.2 | Centered, lockable | 100/1000 |
| 1012.204 | 20 | 0.5 | Centered, lockable | 80/800 |
| 1012.604 | 60 | 1 | Offset, lockable | 60/240 |

SYRINGE TRAYS

nutrisafe2 syringe trays are indicated for the daily preparation of the enteral nutrition for an individual patient. Each tray can contain syringes, lockable caps and sampling device.

Syringes supplied in a transparent tray:



| Code | A-VY syringes | | | Sampling straw Ø 3.0 x 4.2 mm | | Caps packed in a supple blister | Quantity case |
|-----------|---------------|------|--------|-------------------------------|--------|---------------------------------|---------------|
| | Quantity | Size | Type | Quantity | Length | Quantity | |
| 1015.928 | 8 | 20ml | offset | 1 | 20cm | 8 | 20 |
| 1015.966 | 6 | 60ml | offset | 1 | 20cm | 8 | 10 |
| 1015.968 | 8 | 60ml | offset | 1 | 20cm | 8 | 10 |
| 1015.9212 | 12 | 20ml | offset | 1 | 20cm | 16 | 10 |

Syringes supplied with a cardboard display box:



| Code | A-VY syringes | | | Sampling straw Ø 3.0 x 4.2 mm | | Caps packed in a supple blister | Quantity case |
|----------|---------------|------|--------|-------------------------------|--------|---------------------------------|---------------|
| | Quantity | Size | Type | Quantity | Length | Quantity | |
| 1016.928 | 8 | 20ml | offset | 1 | 20cm | 8 | 10 |
| 1016.966 | 6 | 60ml | offset | 1 | 20cm | 8 | 10 |
| 1016.968 | 8 | 60ml | offset | 1 | 20cm | 8 | 10 |

ACCESSORIES

EXTENSION TUBES

STANDARD EXTENSION TUBES

These nutriSafe2 extension tubes keep the feeding tube away from the syringe driver to improve patient's comfort and avoid tube dislodgement.

- Purple-coloured PVC extension tube
- Standard bore for standard milk formula or large bore for fortified milk
- Attached cap to close the tube and avoid contamination.



| Code | Standard extension tube | | | | Quantity |
|---------|-------------------------|-----------|-----------|---------------|----------|
| | Length cm | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | Box/case |
| 368.032 | 30 | 2.5 | 1.5 | 0.62 | 50/300 |
| 368.102 | 100 | 2.5 | 1.5 | 1.86 | 50/300 |
| 368.152 | 150 | 2.5 | 1.5 | 2.74 | 50/300 |
| 369.152 | 150 | 4.0 | 2.8 | 9.4 | 25/300 |
| 369.252 | 250 | 4.0 | 2.8 | 15.55 | 25/150 |

EXTENSION TUBES WITH LATERAL PORT

These nutriSafe2 extension tubes are designed with:

- Funnel proximal end
- Roberts clamp on the proximal section to interrupt the administration of enteral feed while the medication is injected in the lateral port
- Lateral port with male nutriSafe2 connection and attached cap
- Distal end with female nutriSafe2 connection (code 339.0172)
- Distal end with female nutriSafe2 connection and removable Christmas tree adaptor (code 339.0182).



| Code | Extension tube with lateral port | | | | Quantity |
|----------|----------------------------------|-----------|-----------|---------------|----------|
| | Length cm | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | Box/case |
| 339.0172 | 15 | 4.2 | 3.0 | 1.8 | 25/300 |
| 339.0182 | 18 | 4.2 | 3.0 | 2.0 | 25/300 |

EXTENSION TUBES WITH CLAMP

This device gives extended connection to pumps and allows greater freedom of movement for the patients:

- Purple-coloured PVC extension tube
- Roberts clamp
- Attached cap to close the tube.



| Code | Extension tube with clamp | | | | Quantity |
|----------|---------------------------|-----------|-----------|---------------|----------|
| | Length cm | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | Box/case |
| 368.1021 | 100 | 2.5 | 1.5 | 1.86 | 50/300 |
| 368.1521 | 150 | 2.5 | 1.5 | 2.74 | 50/300 |

EXTENSION TUBES FOR GASTROSTOMY BUTTON

Extension tubes for gastrostomy allow to connect a nutriSafe2 syringe or a nutriSafe2 pump administration set to a gastrostomy button for feeding or medication.

Two distal tips are available:

- Straight tip for bolus dose administration
- Curved tip for continuous feeding.



| Code | Extension tubes for gastrostomy button | | | | | Quantity |
|-----------|--|-----------|-----------|---------------|-------------------|----------|
| | Length cm | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | Distal connection | Box/case |
| 5307.1103 | 30 | 6 | 4.2 | 4.5 | Straight | 10/120 |
| 5307.2103 | 30 | 4.2 | 3 | 2.4 | Curved to 90° | 10/120 |
| 5307.1106 | 60 | 6 | 4.2 | 8.7 | Straight | 10/120 |
| 5307.2106 | 60 | 4.2 | 3 | 4.6 | Curved to 90° | 10/120 |

SYRINGE ACCESSORIES

VENTED CAPS

Vented caps are used for nutri**safe**2 syringes (20 and 60ml) to enable gravity feeding and protecting the solution of enteral feed from contamination.



| Code | nutri safe 2 syringe | Syringe capacity | Quantity |
|--------|-----------------------------|-----------------------------|----------|
| | Type | Volume | Box/case |
| 818.22 | C-GON | For 20ml syringes | 50/400 |
| 818.62 | C-GON | For 60ml syringes | 50/400 |
| 818.23 | A-VY | For 20ml syringes | 50/400 |
| 818.63 | A-VY | For 60ml syringes | 50/400 |
| 818.07 | VY-P | For both 20ml/60ml syringes | 30/240 |



MALE CAPS FOR nutri**safe**2 SYRINGES

| Code | Packaging | Quantity (box/case) |
|--------|------------------------|---------------------|
| 828.01 | Individually packed | 100/4800 |
| 828.08 | Pack containing 8 caps | 100/1800 |

These caps for syringe are also available in orange. Please contact your Vygon representative.

HOLDER FOR nutri**safe**2 SYRINGES

Designed to hold 20ml and 60ml nutri**safe**2 syringes in vertical position for gravity feeding.
Easy screw fixation, for attachment to any vertical bar.

Code 367.011


Quantity box/case: 5/120



SAMPLING DEVICES

DRAWING-UP STRAWS

The nutri**safe**2 drawing-up straw is a device connected to a nutri**safe**2 syringe for medication sampling (code 817.052) or to fill the syringe with milk / milk formula (other codes).



| Code | Drawing-up straw | | | | Quantity |
|---------|------------------|-----------|-----------|---------------|----------|
| | Length cm | Int. Ø mm | Ext. Ø mm | Prim. vol. ml | Box/case |
| 817.052 | 5 | 1.2 | 2.0 | 0.11 | 100/1800 |
| 817.152 | 15 | 2.0 | 3.3 | 0.55 | 50/900 |
| 817.153 | 15 | 3.0 | 4.2 | 1.2 | 50/900 |
| 817.154 | 15 | 1.5 | 2.5 | 0.33 | 50/900 |
| 817.202 | 20 | 2.0 | 3.3 | 0.71 | 50/900 |
| 817.203 | 20 | 3.0 | 4.2 | 1.5 | 50/900 |
| 817.302 | 30 | 3.0 | 4.2 | 2.2 | 50/600 |
| 817.502 | 50 | 2.0 | 3.3 | 1.6 | 20/240 |

Drawing-up straws are also available in orange. Please contact your Vygon representative.

DRAWING-UP STRAW WITH FILTER

This drawing-up straw device incorporates:

- A male nutri**safe**2 connector
- A 5µm filter.



| Code | Drawing-up device | | | Quantity |
|---------|-------------------|-----------|-----------|----------|
| | Length cm | Ext. Ø mm | Int. Ø mm | Box/case |
| 817.105 | 10 | 2.5 | 1.5 | 50/900 |

AMPOULE SAMPLING NEEDLE

This needle is designed to aspirate medications from an ampoule to administer the medication orally, via a feeding tube or via a gastrostomy device.

Features:

- Male nutrisafe2 connector
- No filter
- Length: 52 mm
- External diameter: 1.1mm.

Code 817.002
Quantity box/case: 100/1800



LOCKABLE BOTTLE CAP FOR SAMPLING

These nutrisafe2 caps can be used with a large variety of vials to replace the original cap.

Once the cap is in place, the content can only be sampled with a nutrisafe2 syringe for enteral use.

Duration of use up to 30 days.

Non-sterile and individually packed.
Compliant NF H35-49

| Code | Cap | Quantity |
|--------|------|----------|
| | Ø mm | Box/case |
| 820.18 | 18 | 20/1920 |
| 820.20 | 20 | 20/1920 |
| 820.22 | 22 | 20/1920 |
| 820.28 | 28 | 20/1920 |



NEEDLELESS VIAL ACCESS CAP

These devices are designed to perforate vials sealed by a rubber cap.

The contents of the vial can then be aspirated using a nutrisafe2 syringe.

After sampling, the vial can be stored with the nutrisafe2 access cap in place with its attached closing cap. Hence, it is secured against accidental access with a standard syringe or hypodermic needle.

| Code | Compatibility | Quantity |
|---------|---------------|----------|
| | Vial Ø mm | Box/case |
| 821.012 | 14 | 50/400 |
| 821.022 | 20 | 25/200 |



nutrisafe2 spike

Connected to a nutrisafe2 syringe, this device allows medication samplings from vials with a rubber cap. The spike features at the proximal end a male nutrisafe2 connection with cap.

| Code | Spike | Quantity | |
|---------|-----------|---------------|----------|
| | Length mm | Prim. vol. ml | Box/case |
| 281.062 | 31.5 | 0.51 | 10/240 |



OTHER ACCESSORIES

CAP FOR nutrisafe2 FEEDING TUBES

Purple-coloured obturator female cap for nutrisafe2 feeding tubes.
Individually packed.

Code 829.01

Quantity box/case: 100/4800



nutrisafe2 - ENFit ADAPTOR

The nutrisafe2 - ENFit™ adaptor is indicated for enteral use only.
It enables the connection between a female ENFit™ connector (e.g. distal end of an ENFit™ administration pump set) and a male nutrisafe2 connector (e.g. a nutrisafe2 feeding tube).
Priming volume: 0.10ml

Code 368.72

Quantity box/case: 100/1800



ENFit - nutrisafe2 ADAPTOR

The ENFit™ - nutrisafe2 adaptor is indicated for enteral use only.
It enables the connection between a female nutrisafe2 connector (e.g. nutrisafe2 syringe) and a male ENFit™ connector (e.g. a nutriFit feeding tube).
Priming volume : 0.25 ml
Non-sterile

Code 0VEACC510

Quantity units/box: 300/box



griplok® SECUREMENT DEVICE

Designed to secure enteral feeding tubes (and nasal cannulae), this device features :

- Fixation with hydrocolloid adhesive
- griplok® technology that allows re-opening in order to get access to the tube to be cleaned or repositioned.

The griplok® securement device can be cut in half in order to reduce its size and adapt it to the patient's size.

Code 5804.11

Quantity box/case: 50/1000



3-WAY STOPCOCK

The 3-way stopcocks allow delivering medication without disconnecting feed. They incorporate a female nutrisafe2 port to connect to nutrisafe2 feeding tubes and 2 male ports.

Code 5802.02 is chemically more resistant.




| Code | Stopcock | | Prim. vol. ml | Quantity |
|---------|---------------|---------------|------------------|----------|
| | Tested at bar | Material | | Box/case |
| 5802.01 | 10 | Polycarbonate | 0.26 | 50/300 |
| 5802.02 | - | Polyamide | 0.26 | 50/300 |

SYRINGE-SYRINGE TRANSFER DEVICE

Connector to transfer medications contained in a nutrisafe2 syringe into another, without need for intermediary or another container.

Code 368.22

Quantity box/case: 30/2880



FEED SET ADAPTOR

Adaptor to connect pump administration sets to nutrisafe2 feeding tubes.

Code 368.02
Quantity box/case: 30/2880

This adaptor is also available in orange. Please contact your Vygon representative.



TAPER CONNECTOR

Taper connector with nutrisafe2 male connection to administer medications or nutrition into a gastrostomy device.

Code 368.12
Quantity box/case: 30/2880



Y CONNECTOR

The Y connector allows delivering medication without disconnecting feed.

Y connector features nutrisafe2 connections:

- 2 male ports at the proximal end
- 1 female port at the distal end

Priming volume: 0.44ml.

Code 368.42
Quantity box/case: 50/300



ENFit™ / nutriFit®


SAFE ENTERAL FEEDING SYSTEMS
FOR CHILDREN AND ADULTS

FEEDING TUBES

PVC FEEDING TUBES

Adapted for children and adults, the PVC nutri**fit** feeding tubes are intended to be used for short-term enteral feeding (up to 5 days) and feature the following characteristics:

- Made of PVC
- Closed distal tip with 2 lateral eyes (4Fr to 12Fr) and 4 lateral eyes (14Fr to 21Fr)
- Small cap with colour code to quickly identify the size of the tube
- Numerical markings every centimetre:
 - from 5cm to 35cm for codes 364, 361 & 363
 - from 5cm to 70cm for codes 366, 362
- Radiopaque
- Transparent to better visualize the administered solution.



| Code | PVC tube | | | | | Colour code | Quantity |
|---------|-----------|----|-----------|-----------|---------------|-------------|----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | | Box/case |
| 364.047 | 40 | 4 | 1.5 | 0.8 | 0.27 | Purple | 25/300 |
| 364.057 | 40 | 5 | 1.7 | 1.0 | 0.39 | Grey | 25/300 |
| 364.067 | 40 | 6 | 2.0 | 1.2 | 0.54 | Light green | 25/300 |
| 364.087 | 40 | 8 | 2.5 | 1.5 | 0.81 | Blue | 25/300 |
| 361.047 | 50 | 4 | 1.5 | 0.8 | 0.32 | Purple | 50/300 |
| 361.057 | 50 | 5 | 1.7 | 1.0 | 0.47 | Grey | 50/300 |
| 361.067 | 50 | 6 | 2.0 | 1.2 | 0.65 | Light green | 50/300 |
| 361.087 | 50 | 8 | 2.5 | 1.5 | 0.99 | Blue | 50/300 |
| 361.107 | 50 | 10 | 3.3 | 2.0 | 1.7 | Black | 50/300 |
| 363.047 | 75 | 4 | 1.5 | 0.8 | 0.45 | Purple | 25/300 |
| 363.057 | 75 | 5 | 1.7 | 1.0 | 0.67 | Grey | 25/300 |
| 363.067 | 75 | 6 | 2.0 | 1.2 | 0.94 | Light green | 25/300 |
| 363.087 | 75 | 8 | 2.5 | 1.5 | 1.4 | Blue | 25/300 |
| 363.107 | 75 | 10 | 3.3 | 2.0 | 2.5 | Black | 25/300 |
| 363.127 | 75 | 12 | 4.0 | 2.5 | 3.9 | White | 25/300 |
| 366.087 | 90 | 8 | 2.5 | 1.5 | 1.7 | Blue | 25/300 |
| 366.107 | 90 | 10 | 3.3 | 2.0 | 3.0 | Black | 25/300 |
| 366.127 | 90 | 12 | 4.0 | 2.5 | 4.6 | White | 25/150 |
| 366.147 | 90 | 14 | 4.5 | 3.0 | 6.6 | Dark green | 45/180 |
| 366.167 | 90 | 16 | 5.0 | 3.5 | 8.9 | Orange | 45/180 |
| 366.187 | 90 | 18 | 6.0 | 4.0 | 11.7 | Red | 45/180 |
| 362.047 | 125 | 4 | 1.5 | 0.8 | 0.7 | Purple | 25/300 |
| 362.057 | 125 | 5 | 1.7 | 1.0 | 1.06 | Grey | 25/300 |
| 362.067 | 125 | 6 | 2.0 | 1.2 | 1.5 | Light green | 25/300 |
| 362.087 | 125 | 8 | 2.5 | 1.5 | 2.3 | Blue | 25/300 |
| 362.107 | 125 | 10 | 3.3 | 2.0 | 4.0 | Black | 25/300 |
| 362.127 | 125 | 12 | 4.0 | 2.5 | 6.3 | White | 25/150 |
| 362.147 | 125 | 14 | 4.5 | 3.0 | 9.0 | Dark green | 45/180 |
| 362.167 | 125 | 16 | 5.0 | 3.5 | 12.3 | Orange | 45/180 |
| 362.187 | 125 | 18 | 6.0 | 4.0 | 16.0 | Red | 45/180 |

POLYURETHANE FEEDING TUBES

WITHOUT STYLET

Adapted for children and adults, the PUR nutri**fit** feeding tubes are intended to be used for medium-term (1 to 4 weeks) and feature the following characteristics:

- Made of PUR
- Closed distal tip with 2 lateral eyes (4Fr to 12Fr) and 4 lateral eyes (14Fr to 18Fr)
- Small cap with colour code to quickly identify the size
- Numerical markings every centimetre:
 - from 5cm to 35cm for codes 1364, 1361 & 1363
 - from 5cm to 70cm for codes 1366, 1362 & 1365
- Radiopaque line from 4Fr to 6Fr, and ¾ radiopaque tubing from 8Fr
- Transparent to better visualize the administered solution.



| Code | PUR tube without stylet | | | | | Colour code | Quantity |
|----------|-------------------------|----|-----------|-----------|---------------|-------------|----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | | Box/case |
| 1364.047 | 40 | 4 | 1.5 | 0.8 | 0.27 | Purple | 25/300 |
| 1364.057 | 40 | 5 | 1.7 | 1.0 | 0.39 | Grey | 25/300 |
| 1364.067 | 40 | 6 | 2.0 | 1.2 | 0.54 | Light green | 25/300 |
| 1364.087 | 40 | 8 | 2.5 | 1.5 | 0.81 | Blue | 25/300 |
| 1361.047 | 50 | 4 | 1.5 | 0.8 | 0.32 | Purple | 50/300 |
| 1361.057 | 50 | 5 | 1.7 | 1.0 | 0.47 | Grey | 50/300 |
| 1361.067 | 50 | 6 | 2.0 | 1.2 | 0.65 | Light green | 50/300 |
| 1361.087 | 50 | 8 | 2.5 | 1.5 | 0.99 | Blue | 50/300 |
| 1361.107 | 50 | 10 | 3.3 | 2.0 | 1.7 | Black | 50/300 |
| 1363.047 | 75 | 4 | 1.5 | 0.8 | 0.45 | Purple | 25/300 |
| 1363.057 | 75 | 5 | 1.7 | 1.0 | 0.67 | Grey | 25/300 |
| 1363.067 | 75 | 6 | 2.0 | 1.2 | 0.94 | Light green | 25/300 |
| 1363.087 | 75 | 8 | 2.5 | 1.5 | 1.4 | Blue | 25/300 |
| 1363.107 | 75 | 10 | 3.3 | 2.0 | 2.5 | Black | 25/300 |
| 1363.127 | 75 | 12 | 4.0 | 2.5 | 3.9 | White | 25/300 |
| 1366.087 | 90 | 8 | 2.5 | 1.5 | 1.7 | Blue | 25/300 |
| 1366.107 | 90 | 10 | 3.3 | 2.0 | 3.0 | Black | 25/300 |
| 1366.127 | 90 | 12 | 4.0 | 2.5 | 4.6 | White | 25/150 |
| 1362.047 | 125 | 4 | 1.5 | 0.8 | 0.7 | Purple | 25/300 |
| 1362.057 | 125 | 5 | 1.7 | 1.0 | 1.06 | Grey | 25/300 |
| 1362.067 | 125 | 6 | 2.0 | 1.2 | 1.5 | Light green | 25/300 |
| 1362.087 | 125 | 8 | 2.5 | 1.5 | 2.3 | Blue | 25/300 |
| 1362.107 | 125 | 10 | 3.3 | 2.0 | 4.0 | Black | 25/300 |
| 1362.127 | 125 | 12 | 4.0 | 2.5 | 6.3 | White | 25/150 |
| 1362.147 | 125 | 14 | 4.5 | 3.0 | 9.0 | Dark green | 45/180 |
| 1362.167 | 125 | 16 | 5.0 | 3.5 | 12.3 | Orange | 45/180 |
| 1362.187 | 125 | 18 | 6.0 | 4.0 | 16.0 | Red | 45/180 |

WITH STYLET

Adapted for children and adults, the PUR nutriFit feeding tubes are intended to be used for medium-term enteral feeding (1 to 4 weeks) and feature the following characteristics:

- Made of PUR
- Stylet protected by a ball end to facilitate the placement
- Distal tip with 2 lateral eyes (6Fr to 12Fr) and 4 lateral eyes (14Fr to 18Fr)
- Small cap with colour code to quickly identify the size
- Numerical markings every centimetre
 - from 5cm to 35cm for code 1363
 - from 5cm to 70cm for codes 1366 & 1362
- Radiopaque.



| Code | PUR tube with stylet | | | | | | Colour code | Quantity |
|-----------|----------------------|----|-----------|-----------|------------|---------------|-------------|----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Distal tip | Prim. vol. ml | | Box/case |
| 1366.0671 | 90 | 6 | 2.0 | 1.2 | Closed | 1.1 | Light green | 25/150 |
| 1366.0871 | 90 | 8 | 2.5 | 1.7 | Closed | 2.2 | Blue | 25/150 |
| 1366.1071 | 90 | 10 | 3.3 | 2.5 | Closed | 4.6 | Black | 25/150 |
| 1366.1271 | 90 | 12 | 4.0 | 3.2 | Closed | 7.5 | White | 25/150 |
| 1366.1471 | 90 | 14 | 4.5 | 3.5 | Closed | 9.0 | Dark green | 25/150 |
| 1366.1671 | 90 | 16 | 5.3 | 4.0 | Closed | 11.7 | Orange | 20/120 |
| 1366.1871 | 90 | 18 | 6.0 | 4.5 | Closed | 14.8 | Red | 20/120 |
| 1366.0872 | 90 | 8 | 2.5 | 1.7 | Open | 2.2 | Blue | 25/150 |
| 1366.1072 | 90 | 10 | 3.3 | 2.5 | Open | 4.6 | Black | 25/150 |
| 1366.1272 | 90 | 12 | 4.0 | 3.2 | Open | 7.5 | White | 25/150 |
| 1362.0671 | 125 | 6 | 2.0 | 1.2 | Closed | 1.5 | Light green | 25/150 |
| 1362.0871 | 125 | 8 | 2.5 | 1.7 | Closed | 2.9 | Blue | 25/150 |
| 1362.1071 | 125 | 10 | 3.3 | 2.5 | Closed | 6.3 | Black | 25/150 |
| 1362.1271 | 125 | 12 | 4.0 | 3.2 | Closed | 10.3 | White | 25/150 |
| 1362.1471 | 125 | 14 | 4.5 | 3.5 | Closed | 12.3 | Dark green | 25/150 |
| 1362.1671 | 125 | 16 | 5.3 | 4.0 | Closed | 16.0 | Orange | 20/120 |
| 1362.1871 | 125 | 18 | 6.0 | 4.5 | Closed | 20.3 | Red | 20/120 |
| 1362.0872 | 125 | 8 | 2.5 | 1.7 | Open | 2.9 | Blue | 25/150 |
| 1362.1072 | 125 | 10 | 3.3 | 2.5 | Open | 6.3 | Black | 25/150 |
| 1362.1272 | 125 | 12 | 4.0 | 3.2 | Open | 10.3 | White | 25/150 |

SILICONE FEEDING TUBES

WITHOUT WEIGHTED TIP

Adapted for children and adults, the silicone gastro-duodenal nutri**fit** feeding tubes are intended to be used up to 4 weeks and feature the following characteristics:

- Made of Silicone
- Closed distal tip with 2 lateral eyes (4Fr to 10Fr) and 4 lateral eyes (12Fr and 14Fr)
- Small cap with colour code to quickly identify the size
- Numerical markings every centimetre:
 - from 5cm to 35cm for code 2331
 - from 5cm to 70cm for codes 2336 & 2332
- Radiopaque
- Transparent to visualize the administered solution.



| Code | Tube | | | | | Colour code | Quantity |
|----------|-----------|----|-----------|-----------|---------------|-------------|----------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | | Box/case |
| 2331.047 | 50 | 4 | 1.5 | 0.8 | 0.37 | Purple | 10/240 |
| 2331.057 | 50 | 5 | 1.7 | 1.0 | 0.52 | Grey | 10/240 |
| 2331.067 | 50 | 6 | 2.0 | 1.2 | 0.7 | Light green | 10/240 |
| 2331.087 | 50 | 8 | 2.5 | 1.5 | 1.04 | Blue | 10/240 |
| 2336.067 | 90 | 6 | 2.0 | 1.2 | 1.2 | Light green | 10/240 |
| 2332.047 | 125 | 4 | 1.5 | 0.8 | 0.74 | Purple | 10/240 |
| 2332.067 | 125 | 6 | 2.0 | 1.2 | 1.5 | Light green | 10/240 |
| 2332.087 | 125 | 8 | 2.5 | 1.5 | 2.4 | Blue | 10/240 |
| 2332.107 | 125 | 10 | 3.3 | 2.0 | 4.1 | Black | 10/240 |
| 2332.127 | 125 | 12 | 4.0 | 2.5 | 6.2 | White | 10/120 |
| 2332.147 | 125 | 14 | 4.5 | 2.5 | 6.2 | Dark green | 10/120 |

WITH WEIGHTED TIP

Adapted for children and adults, the silicone gastro-duodenal Nutrifit feeding tubes are intended to be used up to 4 weeks and feature the following characteristics:

- Made of Silicone
- Closed distal tip with 2 lateral eyes (6 to 10Fr) and 4 lateral eyes (12Fr & 14Fr)
- Small cap with colour code to quickly identify the size
- Numerical markings every centimetre: from 5cm to 70cm (all sizes)
- Radiopaque
- Transparent to visualize the administered solution.



| Code | Tube | | | | | | Quantity |
|----------|-----------|----|-----------|-----------|---------------|-------------|------------|
| | Length cm | Fr | Ext. Ø mm | Int. Ø mm | Prim. vol. ml | Colour code | Box / case |
| 2395.067 | 125 | 6 | 2.0 | 1.00 | 1.1 ml | Light green | 10 / 240 |
| 2395.087 | 125 | 8 | 2.5 | 1.30 | 1.8 ml | Blue | 10 / 240 |
| 2395.107 | 125 | 10 | 3.3 | 2.00 | 4.1 ml | Black | 10 / 240 |
| 2395.127 | 125 | 12 | 4.0 | 2.50 | 6.2 ml | White | 10 / 120 |
| 2395.147 | 125 | 14 | 4.5 | 2.50 | 6.2 ml | Dark green | 10 / 120 |

SYRINGE TRAYS



Supplied in a blister bag



supplied in a transparent case

| Code | Type of tray | Syringes | Sampling straws | Number of caps | Extension tubes 150cm | Qty/ case |
|------------|------------------|------------|--|----------------|-----------------------|-----------|
| 1025.716 | Transparent case | 6 x 10 ml | 1 Large bore Sampling straw (20cm) | 8 | 6 | 20 |
| 1025.726 | Transparent case | 6 x 20 ml | 2 Sampling straws (20cm): 1 Large & 1 X-Large bore | 8 | 6 | 20 |
| 1015.728 | Transparent case | 8 x 20 ml | 2 Sampling straws (20cm): 1 Large & 1 X-Large bore | 8 | - | 20 |
| 1025.728 | Transparent case | 8 x 20 ml | 2 Sampling straws (20cm): 1 Large & 1 X-Large bore | 8 | 8 | 20 |
| 1015.766 | Transparent case | 6 x 60 ml | 2 Sampling straws (20cm): 1 Large & 1 X-Large bore | 8 | - | 10 |
| 1025.766 | Transparent case | 6 x 60 ml | 2 Sampling straws (20cm): 1 Large & 1 X-Large bore | 8 | 6 | 10 |
| 1015.768 | Transparent case | 8 x 60 ml | 2 Sampling straws (20cm): 1 Large & 1 X-Large bore | 8 | - | 10 |
| 1025.768 | Transparent case | 8 x 60 ml | 2 Sampling straws (20cm): 1 Large & 1 X-Large bore | 8 | 8 | 10 |
| 1015.92127 | Transparent case | 12 x 20 ml | 1 X-Large bore Sampling straw (20cm) | 16 | - | 20 |
| 1015.7281 | Blister bag | 8 x 20 ml | 2 Sampling straws (20cm): 1 Large & 1 X-Large bore | 8 | - | 60 |
| 1015.7681 | Blister bag | 8 x 60 ml | 2 Sampling straws (20cm): 1 Large & 1 X-Large bore | 8 | - | 30 |

ACCESSORIES

EXTENSION TUBES

STANDARD EXTENSION TUBES

These nutrifit extension tubes keep the feeding tube away from the syringe pump or from the feeding pump to improve the patient's comfort.

- Purple-coloured PVC extension tube
- Attached cap to close the tube and avoid contamination depending on the reference.



| Code | Tube | | | | | Quantity |
|---------------------|-----------|-----------|-----------|-------------------|-----|-------------------|
| | Length cm | Ext. Ø mm | Int. Ø mm | Priming volume ml | Cap | Box/case |
| 368.097 | 90 | 2.5 | 1.5 | 1.38 | No | 200/200 |
| 368.107 | 100 | 2.5 | 1.5 | 1.9 | Yes | 45/270 |
| 368.157 or 368.1572 | 150 | 2.5 | 1.5 | 2.8 or 2.3 | No | 45/270 or 200/200 |
| 369.157 | 150 | 4.2 | 3.0 | 9.9 | Yes | 200/200 |

NOTE: Some codes changed, please find below the correlation with the previous codes :
- 368.097 = VZPRO90ENFIT
- 368.1572 = VZPRO150ENFIT
- 369.157 = VZPRO150GDENFIT

EXTENSION TUBES WITH CLAMP

This nutrifit device gives extended connections to pumps and allows greater freedom of movement for ambulatory patients:

- Purple – coloured PVC extension tube
- With Roberts clamp
- Attached cap to close the tube and avoid contamination.

| Code | Tube | | | | Quantity |
|----------|-----------|-----------|-----------|-------------------|----------|
| | Length cm | Ext. Ø mm | Int. Ø mm | Priming volume ml | Box/case |
| 368.1071 | 100 | 2.5 | 1.5 | 1.9 | 45/270 |
| 368.1571 | 150 | 2.5 | 1.5 | 2.8 | 45/270 |



EXTENSION TUBES WITH LATERAL PORT

This nutrifit extension tube is designed with:

- Male ENFit connection at the proximal end
- Roberts clamp on the proximal section to interrupt the administration of enteral feed while the medication is injected in the lateral port
- Lateral port with male ENFit connection and attached cap
- Distal end with female ENFit connection and removable Christmas tree adaptor.

Length (without the Christmas tree adaptor) : 15 cm

Code 0VENSPSSY01
Quantity box : 30



EXTENSION TUBES FOR GASTROSTOMY BUTTON

These nutriFit extension tubes for gastrostomy make it possible to connect a nutriFit syringe or a nutriFit feeding pump set to a gastrostomy button for feeding or medication injection.

Two distal tips are available:

- Straight tip for bolus dose administration
- Curved tip for continuous feeding.
- Non-sterile



| Code | Extension tubes for gastrostomy | | | | | Quantity |
|------------|---------------------------------|--------------|--------------|--------------|----------------------|----------|
| | Length cm | Ext. Ø mm | Int. Ø mm | Lateral port | Distal connection | Box/case |
| 0VENSPMCS3 | 30 | 4.1 | 3 | No | Curved to 90° | 4/40 |
| 0VENSPMCS6 | 65 | 4.1 | 3 | No | Curved to 90° | 4/40 |
| 0VENSPMCY3 | 30 | 4.1 | 3 | Yes | Curved to 90° | 4/40 |
| 0VENSPMCY6 | 65 | 4.1 | 3 | Yes | Curved to 90° | 4/40 |
| 0VENSPMDS3 | 30 | 6.1 | 4.7 | No | Straight | 4/40 |
| 0VENSPMDS6 | 65 | 6.1 | 4.7 | No | Straight | 4/40 |

SYRINGES ACCESSORIES

MALE CAPS FOR nutriFit SYRINGES

These caps are specifically designed for nutriFit syringes. Purple-coloured, it enables instant identification for enteral use.

Compatible with Low Dose Tip syringes.



| Code | Packaging | Quantity box |
|---------|---------------------------|--------------|
| 828.017 | Individually packed | 200 |
| 828.087 | Blister containing 8 caps | 200 |

ENPLUS / ENFit™ CONNECTOR WITH CLAMP

The ENPlus / ENFit™ connector with clamp allows to take a sample of the feeding solution from a female ENPlus nutrition bag to a female ENFit syringe.

Code 0VEACC360
Quantity box: 100



VENTED CAPS

Vented cap are used for nutriFit syringes (20 and 60ml) to enable gravity feeding and protecting the solution of enteral feed from contamination.

Code 818.07
Quantity box: 30/240

NOTE: This code changed, please find here the previous code : VZOPE20/60ENFIT



SAMPLING DEVICES

DRAWING-UP STRAWS

The nutriFit drawing-up straw is a device connected to a nutriFit syringe for medication sampling.



| Code | | Drawing-up straw | | | | | | Quantity | |
|---------|----------|------------------|--------------|--------------|---|------------------|------|----------|---------|
| | | Length cm | Ext. Ø mm | Int. Ø mm | | Prim. vol. ml | | Box/case | |
| 817.057 | 817.0572 | 5 | 2.0 | 1.2 | 1 | 0.13 | 0.10 | 100/1800 | 200/200 |
| 817.157 | | 15 | 3.3 | 2.0 | | 0.57 | | 50/900 | |
| 817.207 | | 20 | 4.1 | 2.8 | | 1.3 | | 200/200 | |
| 817.208 | | 20 | 7.4 | 4.9 | | 3.75 | | 200/200 | |

NOTE: Some codes changed, please find below the correlation with the previous codes :

- 817.207 = VZCAN20ENFIT
- 817.208 = VZCAN20GDENFIT
- 817.0572 = VZCAN5ENFIT

BOTTLE CAPS FOR SAMPLING

These nutriFit caps secure the medicine bottle and allow the sampling with a nutriFit syringe.



| Code | Bottle caps for sampling | | Quantity |
|----------|--------------------------|--|----------|
| | Ø mm | | Box |
| 820.1297 | 12.9 | | 30/240 |
| 820.2047 | 20.4 | | 30/240 |
| 820.2777 | 27.7 | | 30/240 |

NOTE: These codes changed, please find below the correlation with the previous codes :

- 820.1297 = VZBTFSPMENFIT
- 820.2047 = VZBTFSMMENFIT
- 820.2777 = VZBTFSGMENFIT

OTHER ACCESSORIES

CAP FOR FEEDING TUBES

Purple-coloured obturator female cap for nutri**fit** feeding tubes

Code 829.07

Quantity box: 200

NOTE: This code changed, please find here the previous code :VZBOUFENFIT



nutri**safe**2 - ENFit ADAPTOR

The nutri**safe**2 -ENFit™ adaptor is indicated for enteral use only.

It enables the connection between a female ENFit™ connector (e.g. distal end of an ENFit™ administration pump set) and a male nutri**safe**2 connector (e.g. a nutri**safe**2 feeding tube).

Priming volume: 0.10ml

Code 368.72

Quantity box/case: 100/1800



ENFit - nutri**safe**2 ADAPTOR

The ENFit™ - nutri**safe**2 adaptor is indicated for enteral use only.

It enables the connection between a female nutri**safe**2 connector

(e.g nutri**safe**2 syringe) and a male ENFit™ connector (e.g. a nutri**fit** feeding tube).

Priming volume : 0.25 ml

Non-sterile

Code 0VEACC510

Quantity box: 300



3-WAY STOPCOCK

This nutri**fit** 3-way stopcock allows to deliver medication without disconnecting the feed. It incorporates a female nutri**fit** port to connect to nutri**fit** feeding tubes and 2 male ports.

Code 5802.07

Quantity box: 30/240

NOTE: This code changed, please find here the previous code :VZROB3VENFIT



Y CONNECTOR

The nutri**fit** Y connector allows to deliver medication without disconnecting the feed.

Y connector featuring nutri**fit** connections:

- 2 male ports at the proximal end
- 1 female port at the distal end.

Code 368.47

Quantity box/case: 50/300



ENCLEAN BRUSH®

The EnClean Brush® has been specifically designed to fit and aid in the cleaning of the ENFit™ connectors:

- Single-patient use
- Patent protected
- Reusable up to 30 days
- Non-sterile

Code VLENC1252122

Quantity box/case: 10/10



NON-SAFETY
ENTERAL FEEDING
DEVICES


FUNNELED FEEDING TUBES

PVC FEEDING TUBES

WITHOUT WEIGHTED TIP

The funneled PVC feeding tubes are intended to be used for short term use (up to 5 days) and feature the following characteristics:

- Closed or open tip
- Proximal end:
 - Small flexible funnel with cap with colour code for Fr 4 to Fr 10
 - Large flexible funnel without cap for Fr 12 to Fr 21
- Centimetric markings every centimetre and numerical markings every 5cm:
 - from 5cm to 35cm for codes 310 & 311
 - from 5cm to 70cm for code 312
 - from 5cm to 75cm for code 391 & 390
- Radiopaque.
- Transparent to visualize the administered solution.




| Code | Tube | | | | Funnel | | Quantity |
|--------|-----------|----|--------|--------------|--------|-------------|----------|
| | Length cm | Fr | Tip | Lateral eyes | Type | Colour code | Box/case |
| 310.04 | 40 | 4 | Closed | 2 | Small | White | 50/600 |
| 310.05 | 40 | 5 | Closed | 2 | Small | Grey | 50/600 |
| 310.06 | 40 | 6 | Closed | 2 | Small | Green | 50/600 |
| 310.07 | 40 | 7 | Closed | 2 | Small | Yellow | 50/600 |
| 310.08 | 50 | 8 | Closed | 2 | Small | Blue | 50/600 |
| 310.10 | 50 | 10 | Closed | 2 | Small | Black | 50/600 |
| 310.12 | 50 | 12 | Closed | 2 | Large | - | 50/600 |
| 311.04 | 75 | 4 | Closed | 2 | Small | White | 50/600 |
| 311.05 | 75 | 5 | Closed | 2 | Small | Grey | 50/600 |
| 311.06 | 75 | 6 | Closed | 2 | Small | Green | 50/600 |
| 311.08 | 75 | 8 | Closed | 2 | Small | Blue | 50/600 |
| 311.10 | 75 | 10 | Closed | 2 | Small | Black | 50/600 |
| 312.04 | 125 | 4 | Closed | 2 | Small | White | 25/600 |
| 312.06 | 125 | 6 | Closed | 2 | Small | Green | 25/600 |
| 312.07 | 125 | 7 | Closed | 2 | Small | Yellow | 25/600 |
| 312.08 | 125 | 8 | Closed | 2 | Small | Blue | 25/600 |
| 312.09 | 125 | 9 | Closed | 2 | Small | Brown | 25/600 |
| 312.10 | 125 | 10 | Closed | 2 | Small | Black | 25/600 |
| 312.12 | 125 | 12 | Closed | 2 | Large | - | 25/300 |
| 391.06 | 125 | 6 | Closed | 4 | Small | Green | 50/450 |
| 391.08 | 125 | 8 | Closed | 4 | Small | Blue | 50/450 |
| 391.10 | 125 | 10 | Closed | 4 | Small | Black | 50/300 |
| 391.12 | 125 | 12 | Closed | 4 | Large | - | 50/300 |
| 391.14 | 125 | 14 | Closed | 4 | Large | - | 50/300 |
| 391.16 | 125 | 16 | Closed | 4 | Large | - | 50/200 |
| 391.18 | 125 | 18 | Closed | 4 | Large | - | 50/200 |
| 391.21 | 125 | 21 | Closed | 4 | Large | - | 50/200 |
| 390.10 | 100 | 10 | Open | 4 | Small | Black | 50/300 |
| 390.12 | 100 | 12 | Open | 4 | Large | - | 50/300 |
| 390.14 | 100 | 14 | Open | 4 | Large | - | 50/300 |
| 390.16 | 100 | 16 | Open | 4 | Large | - | 50/300 |
| 390.18 | 100 | 18 | Open | 4 | Large | - | 50/200 |
| 390.21 | 100 | 21 | Open | 4 | Large | - | 50/200 |

WITH WEIGHTED TIP

The funneled PVC weighted gastroduodenal feeding tubes feature the following characteristics:

- Weighted metal tip with lateral eyes
- Proximal end:
 - Small flexible funnel with cap with colour code for Fr 4 to Fr 10
 - Large flexible funnel without cap for Fr 12 to Fr 20
- Numerical markings every centimetre from 5cm to 70cm
- Radiopaque
- Transparent to visualize the administered solution.



| Code | Tube | | | Funnel | | Box/case |
|--------|-----------|----|--------------|--------|-------------|----------|
| | Length cm | Fr | Lateral eyes | Type | Colour code | |
| 393.06 | 125 | 6 | 3 | Small | Green | 25/600 |
| 393.08 | 125 | 8 | 3 | Small | Blue | 25/600 |
| 393.10 | 125 | 10 | 3 | Small | Black | 25/600 |
| 392.10 | 125 | 10 | 4 | Small | Black | 50/300 |
| 392.12 | 125 | 12 | 4 | Large | - | 50/300 |
| 392.14 | 125 | 14 | 4 | Large | - | 50/200 |
| 392.16 | 125 | 16 | 4 | Large | - | 50/200 |
| 392.18 | 125 | 18 | 4 | Large | - | 50/200 |
| 392.20 | 125 | 20 | 4 | Large | - | 50/200 |

POLYURETHANE FEEDING TUBES

WITHOUT WEIGHTED TIP

The funneled PUR feeding tubes are intended for medium term use (1 to 4 weeks) and feature the following characteristics:

- Closed tip
- Proximal end:
 - Small flexible funnel with cap with colour code for Fr 4 to Fr 10
 - Large flexible funnel without cap for Fr 12 to Fr 21
- Numerical markings every centimetre:
 - from 5cm to 35cm for codes 1310 & 1311
 - from 5cm to 70cm for code 1391
- Radiopaque
- Transparent to visualize the administered solution.




| Code | Tube | | | Funnel | | Quantity |
|---------|-----------|----|--------------|--------|-------------|----------|
| | Length cm | Fr | Lateral eyes | Type | Colour code | Box/case |
| 1310.04 | 40 | 4 | 2 | Small | White | 25/600 |
| 1310.05 | 40 | 5 | 2 | Small | Grey | 25/600 |
| 1310.06 | 40 | 6 | 2 | Small | Green | 25/600 |
| 1310.08 | 50 | 8 | 2 | Small | Blue | 25/600 |
| 1311.04 | 75 | 4 | 2 | Small | White | 25/600 |
| 1311.05 | 75 | 5 | 2 | Small | Grey | 25/600 |
| 1311.06 | 75 | 6 | 2 | Small | Green | 25/600 |
| 1311.08 | 75 | 8 | 2 | Small | Blue | 25/600 |
| 1391.06 | 125 | 6 | 2 | Small | Green | 50/450 |
| 1391.08 | 125 | 8 | 2 | Small | Blue | 50/450 |
| 1391.10 | 125 | 10 | 4 | Small | Black | 50/300 |
| 1391.12 | 125 | 12 | 4 | Large | - | 50/300 |
| 1391.14 | 125 | 14 | 4 | Large | - | 50/300 |
| 1391.16 | 125 | 16 | 4 | Large | - | 50/200 |
| 1391.18 | 125 | 18 | 4 | Large | - | 50/200 |
| 1391.21 | 125 | 21 | 4 | Large | - | 50/200 |

SILICONE FEEDING TUBES

WITHOUT WEIGHTED TIP

The funneled silicone Levin gastroduodenal feeding tubes feature the following characteristics:

- Proximal end:
 - Small flexible funnel with cap for Fr 6 to Fr 10
 - Large flexible funnel without cap for Fr 12 to Fr 21
- Numerical markings every 5 centimetres:
 - from 20cm to 50cm for Fr 6 to Fr 10
 - from 45cm to 75cm for Fr 12 to Fr 21
- Radiopaque
- Transparent to visualize the administered solution.




| Code | Tube | | | | Funnel | Quantity |
|---------|-----------|----|--------|--------------|--------|----------|
| | Length cm | Fr | Tip | Lateral eyes | Type | Box/case |
| 2391.06 | 125 | 6 | Closed | 4 | Small | 10/240 |
| 2391.08 | 125 | 8 | Closed | 4 | Small | 10/240 |
| 2391.09 | 125 | 9 | Closed | 4 | Small | 10/240 |
| 2391.10 | 125 | 10 | Closed | 4 | Small | 10/240 |
| 2391.12 | 125 | 12 | Closed | 4 | Large | 10/240 |
| 2391.14 | 125 | 14 | Closed | 4 | Large | 10/240 |
| 2391.16 | 125 | 16 | Closed | 4 | Large | 10/240 |
| 2391.18 | 125 | 18 | Closed | 4 | Large | 10/120 |
| 2391.21 | 125 | 21 | Closed | 4 | Large | 10/120 |
| 2316.10 | 125 | 10 | Open | 1 | Small | 10/240 |
| 2316.12 | 125 | 12 | Open | 1 | Large | 10/240 |
| 2316.14 | 125 | 14 | Open | 1 | Large | 10/240 |
| 2316.16 | 125 | 16 | Open | 1 | Large | 10/240 |
| 2316.18 | 125 | 18 | Open | 1 | Large | 10/120 |
| 2316.21 | 125 | 21 | Open | 1 | Large | 10/120 |

WITH WEIGHTED TIP

The funneled silicone weighted gastroduodenal feeding tubes feature the following characteristics:

- Distal end:
 - Brass weighted tip for codes 2394.10 & 2394.12
 - Stainless steel weighted tip for codes 2394.14/16/18/21 & 2395
 - Tungsten weighted tip for code 2396
- Proximal end:
 - Small flexible funnel with cap for Fr 5 to Fr 10
 - Large flexible funnel without cap for Fr 12 to Fr 21
- Numerical markings every 5 centimetres:
 - from 45cm to 75cm for code 2394
 - from 20cm to 50cm for code 2395
 - from 25cm to 75cm for code 2396
- Radiopaque
- Transparent to visualize the administered solution.



| Code | Tube | | | Funnel | Quantity |
|---------|-----------|----|--------------|--------|----------|
| | Length cm | Fr | Lateral eyes | Type | |
| 2396.06 | 90 | 6 | 4 | Small | 25/300 |
| 2396.09 | 125 | 9 | 4 | Small | 25/300 |
| 2396.14 | 125 | 14 | 4 | Large | 25/300 |
| 2396.18 | 125 | 18 | 4 | Large | 25/150 |
| 2395.05 | 125 | 5 | 3 | Small | 25/300 |
| 2395.06 | 125 | 6 | 3 | Small | 25/300 |
| 2395.09 | 125 | 9 | 3 | Small | 25/300 |
| 2394.10 | 125 | 10 | 4 | Small | 10/240 |
| 2394.12 | 125 | 12 | 4 | Large | 10/240 |
| 2394.14 | 125 | 14 | 4 | Large | 10/240 |
| 2394.16 | 125 | 16 | 4 | Large | 10/240 |
| 2394.18 | 125 | 18 | 4 | Large | 10/120 |
| 2394.21 | 125 | 21 | 4 | Large | 10/120 |

JEJUNOSTOMY TRAY

This tray allows the placement of a catheter in the jejunum for patients presenting tumors in the mouth, pharynx, larynx, digestive tract or for patients who have been submitted to surgery in these organs, who cannot be fed by oral nutrition and when the placement of a feeding tube is impossible.



- This jejunostomy tray contains:
- 1 fully radiopaque silicone catheter, length 60 cm \varnothing 2.0 x 3.2 mm
 - 1 Desilet introducer
 - 1 Alene needle Fr 10
 - 1 detachable compression hub (code 2196.30)
 - 1 17G needle – Length 98 mm
 - 2 male Luer-lock/male Luer-lock connectors (code 893.00)

Code: 2300.00
Quantity box/case: 1/20

GASTRIC ASPIRATION FAUCHER GASTRIC TUBES

The Faucher tube features a large diameter tube inserted through the mouth of the patient in the stomach for emergency suctioning and washout. It is used, for example, in the treatment of acute drug intoxication :

- Made of PVC
- Open tip with 4 lateral eyes
- Numerical markings every 5 cm from 15 to 75 cm.



| Code | Length cm | Fr | Funnel | Quantity Box |
|--------|--------------|----|----------|-----------------|
| 351.27 | 150 | 27 | V-shaped | 15 |
| 351.30 | 150 | 30 | V-shaped | 15 |
| 351.33 | 150 | 33 | V-shaped | 15 |
| 351.36 | 150 | 36 | V-shaped | 15 |

DUAL FLOW GASTRIC TUBES

WITHOUT TRAP

The dual flow gastric tubes comprise a principal suction lumen with multiple eyes at the distal end and a secondary ventilation channel with a smaller cross section that communicates with the main suction lumen at the distal tip. A double taper adaptor is provided at the distal end of the main lumen to increase the choice of possible connections. For non-return valve preventing reflux from the lateral port, see page 67

- Made of PVC or PUR
- Closed tip with multiple eyes
- Numerical markings every centimeter
 - from 5cm to 45 cm for 6Fr & 8Fr
 - from 10cm to 65cm for 10Fr
 - from 10cm to 75cm for 12Fr and more
- Radiopaque
- Transparent to visualize the administered solution.



| Code | Tube | | | | Quantity Box/case |
|---------|----------|--------------|----|--------------|----------------------|
| | Material | Length cm | Fr | Lateral eyes | |
| 340.06 | PVC | 60 | 6 | 6 | 10/180 |
| 340.08 | PVC | 60 | 8 | 6 | 10/180 |
| 340.10 | PVC | 90 | 10 | 10 | 10/180 |
| 340.12 | PVC | 120 | 12 | 10 | 10/180 |
| 340.14 | PVC | 120 | 14 | 10 | 10/180 |
| 340.16 | PVC | 120 | 16 | 10 | 10/180 |
| 340.18 | PVC | 120 | 18 | 10 | 10/180 |
| 340.21 | PVC | 120 | 21 | 10 | 10/180 |
| 340.30 | PVC | 120 | 30 | 6 | 15 |
| 340.36 | PVC | 120 | 36 | 6 | 15 |
| 1340.06 | PUR | 60 | 6 | 6 | 10/180 |
| 1340.08 | PUR | 60 | 8 | 6 | 10/180 |
| 1340.10 | PUR | 90 | 10 | 10 | 10/180 |
| 1340.12 | PUR | 120 | 12 | 10 | 10/180 |
| 1340.14 | PUR | 120 | 14 | 10 | 10/180 |
| 1340.16 | PUR | 120 | 16 | 10 | 10/180 |
| 1340.18 | PUR | 120 | 18 | 10 | 10/180 |
| 1340.21 | PUR | 120 | 21 | 10 | 10/120 |

DUAL FLOW GASTRIC TUBES

WITH TRAP: neosump

neosump is a dual flow gastric tube with a collection trap. The tube comprises a main lumen for suction/decompression, while the secondary lumen remains open for ventilation. When a gastric sample is needed, the trap enables collection with accurate measurement of the volume aspirated from the stomach. An additional cap, supplied on the bottom of the trap, allows a safe transport of the sample. A double taper adaptor is provided to increase the choice of suction connections.

- Made of PVC
- Closed tip with multiple eyes
- Numerical markings every centimeter
 - from 5cm to 45 cm for 6Fr & 8Fr
 - from 10cm to 65cm for 10Fr
- Radiopaque
- Transparent to visualize the administered solution.



| Code | Tube | | | Sampling trap | Quantity Box/case |
|---------|--------------|----|--------------|---------------|----------------------|
| | Length cm | Fr | Lateral eyes | | |
| 5382.06 | 60 | 6 | 6 | 35ml | 25/100 |
| 5382.08 | 60 | 8 | 6 | 35ml | 25/100 |
| 5382.10 | 90 | 10 | 10 | 35ml | 25/100 |

The trap is also available separately.
Code : 5534.91
Quantity Box/case : 25/150

ACCESSORIES

STYLET FOR SILICONE FEEDING TUBES

This stylet is designed to make a silicone feeding tube with closed tip slightly more rigid in order to facilitate its insertion:

- in stainless steel covered with a polyurethane sheath coated with hydrogel in order to increase its sliding properties inside the tube
- Proximal end made with an ABS cap in order to facilitate the handling
- Length: 130 cm

Code 2390.101: for 8 Fr to 10 Fr silicone feeding tubes with closed end.
Code 2390.201: for 12 Fr to 21 Fr silicone feeding tubes with closed end.



NON-RETURN VALVE FOR DUAL FLOW GASTRIC TUBE

This non-return valve connects to the lateral port of dual flow gastric tubes (codes: 340/1340) in order to avoid any unintended gastric reflux.



| Code | Material of the body | Material of the valve | Connection | Quantity |
|--------|----------------------|-----------------------|------------|----------|
| | | | | Box/case |
| 852.01 | Polycarbonate | Silicone | Male Luer | 30/1440 |

ENFit™ / nutrifit®

TRANSITION PERIOD ACCESSORIES

⚠ Make sure to use these adaptors for enteral feeding administration only

CHRISTMAS TREE ADAPTORS

A wide range of transition accessories to customize your sets and to help managing the ENFit™ transition period.



| Code | Product | Quantity Box/case |
|-----------|---|-------------------|
| 0VEACC316 | Male Luer-lock / Christmas tree adaptor (with cap, non-sterile) | 50/300 |
| 0VEACC406 | Male ENFit™/ Christmas tree adaptor (without cap, non-sterile) | 50/500 |
| 0VEACC416 | Male ENFit™/ Christmas tree adaptor (with cap, non-sterile) | 50/500 |
| 368.12 | nutri safe 2 / Christmas tree adaptor | 30/2880 |
| 801.00 | Female Luer-lock / Christmas tree adaptor | 30/2880 |

 Make sure to use these adaptors for enteral feeding administration only

FUNNEL ADAPTORS

A wide range of transition accessories to customize your sets and to help managing the ENFit™ transition period.



| Code | Product | Quantity Box/case |
|-----------|---|-------------------|
| 0VEACC500 | Female ENFit™ / funnel adaptor (non-sterile) | 50/500 |
| 0VTACC017 | Female Luer-lock / funnel adaptor (non-sterile) | 300 |
| 800.00 | Male Luer / funnel adaptor | 50/2400 |
| 800.01 | Male Luer-lock / funnel adaptor | 50/2400 |
| 368.02 | Female nutri safe 2 / funnel adaptor | 30/2880 |

 Make sure to use these adaptors for enteral feeding administration only



03/2023 - XXXXXXXX000000X

FOR FURTHER INFORMATION, PLEASE CONTACT: **questions@vygon.com**

The specification shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

VYCON – 5, RUE ADELIN – 95440 ECOUEN – FRANCE

RECEPTION: +33 (0)1 39 92 63 63 – SERVICE CLIENTS FRANCE: +33 (0)1 39 92 63 81

EXPORT CUSTOMER SERVICE: +33 (0)1 39 92 64 15

www.vygon.com

